






A Feedback controller for the RF-KO Method for Slow Extraction for NIMMS Helium Synchrotron

Stamatina Detsi (ECE NTUA)

matina.det@gmail.com

Summer Student Program, Meyrin CERN

Supervisors: R. Taylor, E. Benedetto, M. Vretenar




Extraction

- Slow Extraction
- Fast extraction

In general we push the particles out of their stability triangle

In our case we need constant spill to hit the tumor so we need slow






Slow Extraction

It's the process of extraction the beam little by little in every turn so we create a constant flow of particles from the accelerator (uniform spill).

We want Slow Extraction in order to hit the tumor with a beam that has duration and specific energy.






RFKO Method

In order to achieve slow extraction this we apply a controlled voltage kick.

The kick should be updated every time the extraction rate is different than the expected one.






Our Voltage kick

We have AC, that means we need to take into consideration both Amplitude and Frequency.

So ..


- Amplitude Modulation
- Frequency Modulation





Amplitude Modulation

- Use a PI controller for the Voltage Amplitude.
- If we have more extraction rate \rightarrow decrease
- If we have smaller extraction rate \rightarrow increase






Amplitude Modulation


Use a PI controller for the Voltage Amplitude.

If we have more extraction rate \rightarrow decrease

If we have smaller extraction rate \rightarrow increase




PID Controller






$e(t)$: error function



PID Controller



Our Voltage kick



Number of extracted particles every 10.000 turns





Frequency Modulation

To extract the particles we need to put hit the in their tune frequency


We don't have only one tune so we need a way to affect all particles, so we change the voltage constantly, with:

1. Random numbers
2. Sinwave in resonance
3. Sinwave in resonance with same random noise
4. Modulator signal
5. Chirp Method
6. BPSK

To extract the particles we need to put hit the in their tune frequency. We dont have only one tune so we need a way to effect all particles, so we change the voltage all the time, with:


Random

Commonly used as white-noise




Sine wave 1:

At frequency Q_x of beam




Sine wave 2:

Random noise matching beam tune distribution




Modulator signal



Chirp Method

Used in PS and HIMAC





My Simulation

PI Controller

My error

Step Parameter

Calculate Kick

Frequency Modulation



Tracking

Xsuite

Helium Synchrotron

Blackbox

Wrote a Beam Interact function to calculate number of particles extracted per window (defined by the electrostatic septum), and fed this into your controller

Controller gives a increase to the px coordinate of the 6D particle beam




First Approach: Sine Wave 1

Manually changing the parameters (K_p , K_i , setpoint) of PI controller


-> found we have extraction only while our voltage is increasing

-> we have always ONLY 60-80 particles extracted until amplitude too big


-> it looks the same either for 10^3 or 10^5 turns



Modulator Signal example - given to the controller



Modulator Signal - Controller Result



Realpoint-0.1*Setpoint



Fine tuning

Chirp Method


Manual changing the parameters of pi controller

-> all the particles are extracted

-> it looks like the controller is trying to reach the setpoint


-> we have big fluctuation maybe because the parameters are not optimized

Chirp Signal example - given to the controller






Chirp Signal example - given to the controller






Chirp Signal example - given to the controller



Chirp Signal controller result




Realpoint-Setpoint ~ 13%

The importance of Frequency Modulation:

- High amplitude = beam blow-up




Initial Beam




Septum at -55 mm
'extracts' particles
above this amplitude

Final Beam Case 1
Kick too high



Final Beam Case 2
Extracted beam as
expected




The importance of Frequency Modulation:


- High amplitude = beam blow-up




Initial Beam



Final Beam Case 1
Kick too high



Final Beam Case 2
Extracted beam as expected




Septum at -55 mm
'extracts' particles
above this amplitude



RESULTS


Chirp Result - 2000 particles



Question: Does it scale with number of particles?

NIMMS needs 10^{10} particles


Chirp Result - 10000 particles



Need to
re-optimize!

5x more
particles, so
need 5x faster
rate (setpoint)

Chirp Result - 10000 particles - 5x setpoint



Reduced error
and oscillation
size


Spill Characteristic: Probability Mass Function




We are trying to control how many particles are extracted per window of 500 turns. This effect is related to a Poisson distribution



Chirp Result - 2000 particles



Chirp Result - 10000 particles



Spill Characteristic: Duty Factor

Why is it important? Only for
comparing


It means nothing!

Measured $F_{av} = 0.88758$

$$\frac{\sum (N_p/T)^2}{\sum ((N_p/T)^2)}$$




Duty Factor



Attention:

No standardized definition of duty factor!
Direct comparisons are limited!

Chirp Result - 10000 particles + 5x setpoint



Measured $F_{av} = 0.88758$

Conclusions:



Now on:

- Rewriting it in C in Xsuite rather than Python to be faster
- This tools will be used to instantly extract from every Helium Synchrotron setting (num. particles, emittance, sextupole strength, tune distance etc.)
- Before, each RF-KO was hand-optimised from a 8-variable function
- Already have other interested users from Xsuite community to make it into a full element

Further improvement can be made:

- D part of PID controller
- Feedforward System
- Further optimisation of parameters





Extra Slides

Helium Synchrotron Extraction Details

Simulations performed in Xsuite. Whole accelerator is a 'blackbox', but uses the following specifications:

- Resonant sextupole & injection septa in straight-section 1
- Magnetic extraction septa in straight-section 2
- Electrostatic septum in straight-section 3

For simplicity, simulation tracking starts and ends at electrostatic septum, using stable beam that has already circulated for 10,000 turns.



Tunes (x,y)	2.655, 0.705
Chroma (x,y)	-6.382, -3.273
Emittance_n (x,y)	5, 1 [mm.mrad]
Res sextupole k2	-1
ES Septum X_lim	-55 [mm]
delta p/p	0

